

Medis® Suite MR

Hardware Requirements

Recommended Systems			
	Standalone	Client	Server
Operating System	Windows 7, 8.1, 10 64 bit	Windows 7, 8.1, 10 64 bit	Windows 7, 8.1, 10 Windows 2012 (R2), Windows 2016 64 bit
Processor	Octa Core CPU	Octa Core CPU	Octa Core CPU
Memory	16GB+	16GB+	16GB+
Hard Drive	Hard Drive 1TB SSD	Hard Drive 250GB+ SSD	Hard Drive 1: 250GB+ SSD Hard Drive 2: Data Storage: 1TB*
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16 NVIDIA CUDA compatible	Dedicated OpenGL compatible 512MB PCIe x16 NVIDIA CUDA compatible	Dedicated OpenGL compatible 512MB PCIe x16 NVIDIA CUDA compatible
Display	22" widescreen 1920x1080 Resolution	22" widescreen 1920x1080 Resolution	-
<p>Network connection \geq 1000 Mbit/s and latency \leq 1 msec recommended for client server connections.</p> <p>* Recommended for server up to 5 clients and 200 – 500 concurrently available studies. Size hard drive depends on the size and amount of studies.</p>			

Minimum Systems

	Standalone	Client	Server
Operating System	Windows 7, 8.1, 10 64 bit	Windows 7, 8.1, 10 64 bit	Windows 7, 8.1, 10 Windows 2012 (R2), Windows 2016 64 bit
Processor	Quad Core CPU	Quad Core CPU	Quad Core CPU
Memory	8GB	8GB	8GB
Hard Drive	Hard Drive: 250GB+*	Hard Drive: 250GB+*	Hard Drive: 250GB+*
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16
Display	22" widescreen 1920x1080 Resolution	22" widescreen 1920x1080 Resolution	-

Network connection \geq 100 Mbit/s and latency \leq 1 msec recommended for client server connections.

* Hard drive depends on the size and amount of studies, in case of spinning disk: 7200 rpm or more.

Medis Medical Imaging Systems bv

Schuttersveld 9, 2316 XG Leiden
P.O. Box 384, 2300 AJ Leiden, The
Netherlands

P +31 71 522 32 44

F +31 71 521 56 17

E sales@medis.nl

Medis Medical Imaging Systems, Inc.

9360 Falls of Neuse Road, Suite 103
Raleigh, NC 27615-2484, USA

P +01 (919) 278 7888

F +01 (919) 847 8817

E us-sales@medis.nl

© 2019, Medis Medical Imaging Systems bv
8.21.230.32.1