

# Medis® Suite CT

## Hardware Requirements

Recommended Systems			
	Standalone	Client	Server
Operating System	Windows 10 64 bit	Windows 10 64 bit	Windows 10 Windows server 2012 (R2), Windows server 2016 64 bit, Windows server 2019 64 bit
Processor	Quad Core CPU	Quad Core CPU	Quad Core CPU
Memory	8GB+	8GB+	8GB+
Hard Drive	Hard Drive 500GB+ SSD*,**	Hard Drive 250GB+ SSD **	Hard Drive 1: System 250GB+ SSD*,**
	-	-	Hard Drive 2: Study data: 500GB+ *,**
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16 **
Display	22“ widescreen 1920x1080 Resolution	22“ widescreen 1920x1080 Resolution	-
Mouse	3-button mouse***	3-button mouse***	3-button mouse***

Network connection >= 100 Mbit/s and latency <= 1 msec recommended for client server connections  
 \* Hard drive depends on the size and amount of studies  
 \*\* When the software GUI is not used (pure servers), no dedicated graphics card is required.  
 \*\*\* Separate scroll wheel and middle mouse button is recommended.

## Minimum System

	<b>Standalone</b>	<b>Client</b>	<b>Server</b>
Operating System	Windows 10 64 bit	Windows 10 64 bit	Windows 10 Windows server 2012 (R2), Windows server 2016 64 bit Windows server 2019 64 bit
Processor	Quad Core CPU	Quad Core CPU	Quad Core CPU
Memory	4GB	4GB	4GB
Hard Drive	Hard Drive 250GB+ *	Hard Drive 250GB+ *	Hard Drive 250GB+ *
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16 **
Display	22" widescreen 1920x1080 Resolution	22" widescreen 1920x1080 Resolution	-
Mouse	3-button mouse***	3-button mouse***	3-button mouse***

Network connection  $\geq 100$  Mbit/s and latency  $\leq 1$  msec recommended for client server connections

\* Hard drive depends on the size and amount of studies

\*\* When the software GUI is not used (pure servers), no dedicated graphics card is required.

\*\*\* Separate scroll wheel and middle mouse button is recommended.

### Medis medical imaging systems bv

Schuttersveld 9 2316 XG Leiden  
 P.O. Box 384 2300 AJ Leiden, the Netherlands  
 Tel. +31 71 522 32 44  
 Fax + 31 71 521 56 17  
[sales@medisimaging.com](mailto:sales@medisimaging.com)  
[www.medisimaging.com](http://www.medisimaging.com)

### Medis medical imaging systems, Inc.

9360 Falls of Neuse Road, Suite 103  
 Raleigh, NC 27615-2484, USA  
 Tel. +01 (919) 278 7888  
 Fax +01 (919) 847 8817  
[us-sales@medisimaging.com](mailto:us-sales@medisimaging.com)