

Minimum System

	Standalone	Client	Server
Operating System	Windows 10 64 bit	Windows 10 64 bit	Windows 10 Windows MUMU 2012 (R2), Windows MUMU 2016 64 bit Windows server 2019 64 bit
Processor	Quad Core CPU	Quad Core CPU	Quad Core CPU
Memory	4GB	4GB	4GB
Hard Drive	Hard Drive 250GB+ *	Hard Drive 250GB+ *	Hard Drive 250GB+ *
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16 **
Display	22" widescreen 1920x1080 Resolution	22" widescreen 1920x1080 Resolution	-
Mouse	3-button mouse***	3-button mouse***	3-button mouse***

Network connection \geq 100 Mbit/s and latency \leq 1 msec recommended for client server connections

* Hard drive depends on the size and amount of studies

** When the software GUI is not used (pure servers), no dedicated graphics card is required.

*** Separate scroll wheel and middle mouse button is recommended.

Medis medical imaging systems bv
Schuttersveld 9 2316 XG Leiden
P.O. Box 384 2300 AJ Leiden, the Netherlands
Tel. +31 71 522 32 44
Fax + 31 71 521 56 17
sales@medisimaging.com
www.medisimaging.com

Medis medical imaging systems, Inc.
9360 Falls of Neuse Road, Suite 103
Raleigh, NC 27615-2484, USA
Tel. +01 (919) 278 7888
Fax +01 (919) 847 8817
us-sales@medisimaging.com