

Medis® Suite CT Hardware Requirements

Recommended Systems			
	Standalone	Client	Server
Operating System	Windows 11 64 bit	Windows 11 64 bit	Windows servers 2019, 2022 64 bit
Processor	Quad Core CPU	Quad Core CPU	Quad Core CPU
Memory	8GB+	8GB+	8GB+
Hard Drive	500GB+ SSD*	250GB+ SSD*	System: 250GB+ SSD* Study data:500GB+*
Free disk space	10GB**	10GB**	10GB**
Graphics Card	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16	Dedicated OpenGL compatible 512MB PCIe x16 ***
Display	22" widescreen 1920x1080 Resolution	22" widescreen 1920x1080 Resolution	-

Network connection ≥ 100 Mbit/s and latency ≤ 1 msec recommended for client server connections

* Hard drive depends on the size and amount of studies

** Should be available after Medis Suite is installed, and does not take disk space for repository data into account.

*** When the software GUI is not used (pure servers), no dedicated graphics card is required.

Minimum System

	Standalone	Client	Server
Operating System	Windows 10 64 bit	Windows 10 64 bit	Windows servers 2012 (R2), 2016
Memory	8GB	8GB	8GB
Hard Drive	Hard Drive 250GB+ *	Hard Drive 250GB+ *	Hard Drive 250GB+ *
Memory	10GB**	10GB**	10GB**
Display	1.3 Megapixel resolution	1.3 Megapixel resolution	-

Network connection \geq 100 Mbit/s and latency \leq 1 msec recommended for client server connections

* Hard drive depends on the size and amount of studies

** Should be available after Medis Suite is installed, and does not take disk space for repository data into account.

Medis Medical Imaging Systems bv
Schuttersveld 9 2316 XG Leiden
P.O. Box 384 2300 AJ Leiden, the
Netherlands Tel. +31 71 522 32 44
Fax + 31 71 521 56 17
sales@medisimaging.com
www.medisimaging.com

Medis Medical Imaging Systems, Inc.
9360 Falls of Neuse Road, Suite 103
Raleigh, NC 27615-2484, USA
Tel. +01 (919) 278 7888
Fax +01 (919) 847 8817
us-sales@medisimaging.com