

# Medis Suite XA

## All-in-one solution for anatomical & functional cardiovascular X-ray quantitative analyses

Medis Suite XA is a cardiovascular post-processing X-ray imaging software that offers 2D and 3D analysis of the coronary arteries, peripheral arteries and of the ventricles. 7 modules are available in Medis Suite XA, different module combinations can be selected based on your clinical needs.

### FOR QUANTIFIED CLINICAL DECISIONS

- Quantifies lesion anatomy for coronary and peripheral arteries
- Offers measurement of heart function

### WORKFLOW EFFICIENCY

- One simple workflow for all applications, easy acquisition and analysis

### ROBUST PLATFORM

- Robust Medis Suite platform, which integrates browsing, analysing and reporting

### WHICH MODULES ARE PRESENT IN MEDIS SUITE XA ?

- **ANATOMY**
  - 3D QCA Single Vessel
  - 3D QCA Bifurcation
  - 2D QCA Single Vessel
  - 2D QVA Single Vessel
- **HEART FUNCTION**
  - LV Monoplane
  - LV Biplane
  - RV Biplane

### WHAT IS NEW IN MEDIS SUITE XA?

#### QCA

- Fast and auto selection of best ED frame based on AI
- Reference diameter now available in the Lesion tables
- AI based automatic matching on two views
- Additional stent results in 2D QCA


#### LV biplane

- New regression formula now available for adult patients



**AI DEEP LEARNING  
IN MEDIS SUITE XA**

## ALL-IN-ONE SOLUTION FOR YOUR DAILY XA ROUTINE



## FOR QUANTIFIED CLINICAL DECISIONS

- Validated results
- Applicable for coronary bifurcation and peripheral vessels
- For analysis of heart function


## WORKFLOW EFFICIENCY

- One streamlined workflow for all modules
- Only standard angiograms needed
- Independent of X-ray Vendor


## ROBUST PLATFORM

- Medis Suite platform offers role-based access control
- Integrated viewing, analyzing and reporting
- Offers repository management and data handling


## MEDIS SUITE XA OFFERS




2D QCA Single Vessel




3D QCA Single Vessel




LV Monoplane



3D QCA Bifurcation



2D QCA Single Vessel



LV/RV Biplane

### Legal Statements

Medis Suite and QAngio XA have market clearance for USA, Canada, EU, Australia, Japan and South Korea.

QAngio XA 3D (incl. QFR) has market clearance for the USA, Canada, EU, Australia, Brazil, Malaysia, Singapore, Indonesia and South Korea. Market clearance in Japan is pending.

### Medis Medical Imaging Systems BV

Schuttersveld 9, 2316 XG Leiden

P.O. Box 384, 2300 AJ Leiden, The Netherlands

P +31 71 522 32 44 F +31 71 521 56 17 E sales@medisimaging.com

### Medis Medical Imaging Systems Inc.

9360 Falls of Neuse Road, Suite 103

Raleigh, NC 27615-2484, USA

P +01 (919) 278 7888 F +01 (919) 847 8817 E us-sales@medisimaging.com

© 2023, Medis

Medical Imaging

Systems BV

8.22.200.40.6