

Medis Suite MR

Hardware Requirements

April 7, 2025

Recommended System

	Standalone	Client	Server
Operating System	Windows 11 64 bit	Windows 11 64 bit	Windows Server 2019, 2022 64 bit
CPU	8 core	8 core	8 core
Memory	16 GB+	16 GB+	16 GB+
Hard drive	1 TB SSD*	250 GB+ SSD*	System: 250 GB+ SSD* Data Storage: 1TB*
Network	Speed >= 100 Mbps Latency <= 1 ms	Speed >= 100 Mbps Latency <= 50 ms	
Free disk space	10 GB**	10 GB**	10 GB**
Display	22" wide screen 1920 x 1080 px resolution	22" wide screen 1920 x 1080 px resolution	22" wide screen 1920 x 1080 px resolution

Recommended for server up to 5 clients and 200 – 500 concurrently available studies.

* Hard drive depends on the size and amount of studies. In case of spinning disk: 7200 rpm or more.

** Should be available after Medis Suite is installed, and does not take disk space for repository data into account.

Medis Suite MR

Hardware Requirements

April 7, 2025

Minimum System

	Standalone	Client	Server
Operating System	Windows 10 64 bit	Windows 10 64 bit	Windows Server 2012 (R2), 2016
CPU	4 core	4 core	4 core
Memory	8 GB	8 GB	8 GB
Hard drive	Hard Drive 250GB+*	Hard Drive 250GB+*	Hard Drive 250GB+*
Network	-	Speed >= 100 Mbps Latency <= 1 ms	Speed >= 100 Mbps Latency <= 1 ms
Free disk space	10 GB**	10 GB**	10 GB**
Display	1280 x 1024 (4:3) or 1600 x 900 (16:9)	1280 x 1024 (4:3) or 1600 x 900 (16:9)	-

Recommended for server up to 5 clients and 200 – 500 concurrently available studies.

* Hard drive depends on the size and amount of studies. In case of spinning disk: 7200 rpm or more.

** Should be available after Medis Suite is installed, and does not take disk space for repository data into account.